

Urology

A Pediatric Focus



220015, г. Минск, ул. Пономаренко 35а, оф.037
+ 375 17 207 13 67 (/) +375 33 345-18-59 (MTS)
+ 375 29 593 48 56 (Velcom) e- mail: info@medte.by
medte.by

COOK®



COOK® is one of the world's best-known and most respected names in the medical device industry. Closely allied with top medical research facilities, teaching hospitals and leading specialists around the world, Cook continues to provide innovative, physician-driven, quality medical devices to health care providers who will accept nothing but the best for their patients.

Our pediatric urology products continue to provide cutting-edge technology to practitioners who demand better outcomes with fewer complications.


Your high standards in pediatric patient care are enabled by the integrity of Cook. We have provided this assurance for over 40 years and will continue to do so into the future.

We hope you'll be pleased...



Please refer to our Focus Catalogue to view our range of Urology products

CONTENTS	PAGE
PEDIATRIC URETERAL CATHETERS	1
WHISTLE TIP URETERAL CATHETERS	1
ANGLED TIP URETERAL CATHETER SET	2
FLEXOR & FLEXOR DL URETERAL ACCESS SHEATH	2
DOUBLE PIGTAIL PEDIATRIC STENT SETS	3
TARKINGTON URETHRAL STENT SET	3
ZAONTZ URETHRAL STENTS	4
PERCUTANEOUS PIGTAIL NEPHROSTOMY SETS	4
KOYLE DIAPER STENTS	5
PEDIATRIC NEPHROSTOMY CATHETER/STENT SETS	5
SALLE INTRAOPERATIVE PYELOPLASTY STENT SETS	6
MINI-PERC™ ENTRY SET	7
SILICONE FOLEY	8
PEDIATRIC URINARY DRAIN BAG SET	8



PEDIATRIC URETERAL CATHETERS

Used for drainage and retrograde pyelogram in pediatric patients. A stiffening stylet is included to aid in catheter placement. Supplied sterile in peel-open packages. Intended for one-time use.




STIFFENING STYLET
.010 inch stainless steel



CATHETER
TFE

CATHETER

ORDER NUMBER	French Size	Length
022230	3.0	70 cm



WHISTLE TIP URETERAL CATHETERS

Used for drainage and retrograde pyelogram. Supplied sterile in peel-open packages. Intended for one-time use.



CATHETER
Radiopaque polyvinylchloride




ORDER NUMBER	French Size	Length
022103	3.0	70 cm
022104	4.0	70 cm



ANGLED TIP URETERAL CATHETER SET


Used for directing a flexible tipped wire guide in the ureter. The angled tip allows passage of the catheter past a ureteral stone in preparation for E.S.W.L. Supplied sterile in peel-open packages. Intended for one-time use.



ORDER NUMBER	French Size	Length	Remarks
023104	4.0	70 cm	.028 inch (0.71 mm) TFE Wire guide included, tapered end

FLEXOR & FLEXOR DUAL LUMEN URETERAL ACCESS SHEATHS

Used for ureteral dilation and to provide a continuous working channel for the introduction of endoscopes and instruments during ureteral access procedures. Once in position, the Flexor sheath protects the ureter during repeated instrument exchanges, reducing the potential for trauma. The continuous channel also protects delicate instruments and scopes from damage, minimizing costly repairs.



FLEXOR URETERAL ACCESS SHEATH

Illustrated Above

ORDER EXAMPLE: FUS-120013 (12 Internal Diameter and 13 cm length)

ORDER NUMBER	Internal Diameter	Length
FUS-095013	9.5	13 cm
FUS-120013	12	13 cm




FLEXOR DUAL LUMEN URETERAL ACCESS SHEATH

Not Illustrated

ORDER EXAMPLE: FUS-120013-DL (12 Internal Diameter and 13 cm length)


ORDER NUMBER	Internal Diameter	Length
FUS-095013-DL	9.5	13 cm
FUS-120013-DL	12	13 cm





DOUBLE PIGTAIL PEDIATRIC STENT SETS

Used for temporary internal drainage from the ureteropelvic junction to the bladder in pediatric patients. Supplied sterile in peel-open packages. Intended for one-time use. **CAUTION: Periodic evaluation is advised; polyurethane, C-Flex®, and Sof-Flex® stents must not remain indwelling more than six (6) months.**



ORDER NUMBER	French Size	Length	Remarks
0393(08,10,12,14)	3.0	08,10,12,14 cm	Sof-Flex Stents not available with sideports or tether (as illustrated)

ORDER EXAMPLE: 039310 = 3.0 Fr. Stent, 10 cm long, .018 inch Cope Mandril Wire Guide, 4.0 Fr. Stent positioner

0363(08,10,12,14)	3.7	08,10,12,14 cm	C-Flex Set includes .028 inch (0.71 mm) PTFE coated wire guide, Ureteral Catheter (4.0 Fr., 70 cm long and positioner (3.7 Fr., 50 cm long))
-------------------	-----	----------------	--

ORDER EXAMPLE: 036312 = 3.7 Fr. Stent, 12 cm long, .028 inch straight Safety Wire Guide, 3.7 Fr. Stent positioner




TARKINGTON URETHRAL STENT SET

Used for urethral stenting and to provide postoperative drainage following hypospadias repair. The stent is designed to be stabilized in the urethra without use of sutures at the meatus. Supplied sterile in peel-open packages. Intended for one-time use.




ORDER NUMBER	French Size
032408	8.0



ZAONTZ URETHRAL STENTS

Used for stenting the urethra after hypospadias or epispadias repair and to provide post-operative drainage of the bladder. Supplied sterile in peel-open packages. Intended for one-time use.




STENT
Contains C-Flex® TPE

ORDER NUMBER	French Size	Length
032306	6.0	12 cm
032308	8.0	12 cm
032310	10.0	12 cm



PERCUTANEOUS PIGTAIL NEPHROSTOMY SETS

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage. Supplied sterile in peel-open packages. Intended for one-time use. **CAUTION: Periodic evaluation is advised; catheter must not remain indwelling more than four (4) months.**



ORDER NUMBER	French Size	Length	Remarks
080006	6.0	15 cm	Set consists of items shown above



KOYLE DIAPER STENTS

Used for stenting the urethra after hypospadias or epispadias repair and to provide postoperative drainage of the bladder. The externalized portion of the stent is of sufficient length to allow drainage into a second or outside diaper, keeping the inside diaper dry during the healing process. Supplied sterile in peel-open packages. Intended for one-time use.




ORDER NUMBER	French Size	Length
KDSCF-060050	6.0	50 cm
KDSCF-080050	8.0	50 cm



PEDIATRIC NEPHROSTOMY CATHETER/STENT SETS

Used to prevent postoperative stricture and to allow renal drainage after an endopyelotomy. Supplied sterile in peel-open packages. Intended for one-time use. **CAUTION: Periodic evaluation is advised; catheter must not remain indwelling more than four (4) months.**



Silicone




ORDER NUMBER	French Size	Stent Segment Length	Accepts Wire Guide Diameter
080208-56	8.2	5 cm	.038 inch (0.97mm)
080208-57	8.2	7 cm	.038 inch (0.97mm)



SALLE INTRAOPERATIVE PYELOPLASTY STENT SETS

Used to establish and maintain nephrostomy drainage following pyeloplasty. Supplied sterile in peel-open packages. Intended for one-time use. **CAUTION: Periodic evaluation is advised; stent must not remain indwelling more than four (4) weeks.**



ORDER NUMBER	French Size	STENT Usable Length	CATHETER Length	POSITIONING STYLET Length	WIRE GUIDE	
					Diameter	Length
SIPSF-040018-59	4.0	12-18 cm	59 cm	12 cm	.025 inch (0.64 mm)	145 cm
SIPSF-050018-59	5.0	12-18 cm	59 cm	12 cm	.038 inch (0.97 mm)	145 cm



MINI-PERC™ ENTRY SET

Used to gain percutaneous access into the renal pelvis to treat stone disease. The dual lumen catheter and the Peel-Away® introducer are sized for use in pediatric patients. Supplied sterile in peel-open packages. Intended for one-time use.



BENTSON WIRE GUIDE

.035 inch (0.89 mm) diameter PTFE coated stainless steel 145 cm long with 15 cm flexible tip




AMPLATZ ULTRA STIFF WIRE GUIDE

.035 inch (0.89 mm) diameter PTFE coated stainless steel 145 cm long with 8 cm flexible tip



FASCIAL INCISING NEEDLE
18 gage stainless steel 5 cm long




PEEL-AWAY® SHEATH
13.0 French radiopaque TFE 11.5 cm long

INTRODUCER
11.0 French radiopaque polyethylene 13 cm long



DUAL LUMEN ACCESS CATHETER
10.0 French polyurethane 40 cm long

ORDER NUMBER	Remarks
080000-DOCIMO	SET consists of items shown above.



SILICONE FOLEY

Used to provide drainage of urine from the bladder. All-silicone construction provides maximum softness and biocompatibility for enhanced patient comfort. Supplied sterile in peel-open packages. Intended for one-time use.




ORDER NUMBER	French Size	Length	Inflation Volume	Remarks
028506-OE	6.0	23.5 cm	1 ml	Open end; for neonates and non-ambulatory pediatrics.
028506-OE-21	6.0	23.5 cm	1 ml	Open end; with .021 inch diameter SlipCoat™ wire guide. For neonates and non-ambulatory pediatrics.
028508-OE	8.0	23.5 cm	3 ml	Open end



PEDIATRIC URINARY DRAIN BAG SET

Used to provide pediatric patients an appropriately sized disposable system for urine drainage and collection. Supplied non-sterile. Intended for one-time use.



ORDER NUMBER	Remarks
UDPS-200	SET consists of items shown above.

COOK®

Ритейлер в Беларуси - ООО "Медицинские технологии и сервис"
220015, г. Минск, ул. Пономаренко 35а, оф.037
+375 17 207 13 67 (Тел/факс) +375 33 345-18-59 (MTS)
+375 29 593 48 56 (Velcom) e-mail: info@medte.by
medte.by